

HMK15 Humidity Calibrator
With Two Standard Chambers

Vaisala Humidity Calibrator		HMK15											PRICE
Salt chamber	no additional salt chamber	A											
	<i>an additional salt chamber (19766HM)</i>	<i>B</i>											
LiCl salt	no LiCl	0											
	<i>LiCl (11 %RH) for one chamber (19729HM)</i>	<i>1</i>											
	<i>LiCl (11 %RH) for two chambers</i>	<i>2</i>											
MgCl2 salt	no MgCl2	0											
	<i>MgCl2 (33%RH) for one chamber (19730HM)</i>	<i>1</i>											
	<i>MgCl2 (33%RH) for two chambers</i>	<i>2</i>											
NaCl salt	no NaCl	0											
	<i>NaCl (75 %RH) for one chamber (19731HM)</i>	<i>1</i>											
	<i>NaCl (75%RH) for two chambers</i>	<i>2</i>											
K2SO4 salt	no K2SO4	0											
	<i>K2SO4 (97 %RH) for one chamber (19732HM)</i>	<i>1</i>											
	<i>K2SO4 (97 %RH) for two chambers</i>	<i>2</i>											
Ion exchanged water	no water	A											
	<i>water for three chambers (19767HM)</i>	<i>B</i>											
Carrying bag	no bag	0											
	<i>carrying bag (HM27033)</i>	<i>1</i>											
Thermometer	with mercury (19728HM)	1											
	with red capillary liquid (25130HM)	2											
Manual	no manual	A											
	English	B											
	French	C											
										TOTAL			
										QTY			
										TOTAL VALUE			

Selections in bold are included in the prices of the basic versions.
Selections in italic are available at an extra price.
 Shaded sections indicate typical selections.

Example of order code with typical settings:

HMK15	A	1	0	1	0	A	0	1	B
--------------	----------	----------	----------	----------	----------	----------	----------	----------	----------

Bulk Salt Accessories:	Description	Order Code	Size	PRICE
	Lithium Chloride	5236	250 g	
	Sodium Chloride	4102	1 kg	
	Potassium Sulfate	2827	1 kg	
Additional Accessories:	Salt Chamber	19766HM		